

Future suitable disinfectants

Future	Disinfectant	Manufacturer /Supplier
✓	Meiseptol foam pure	B. Braun
✓	Bacillol AF	BODE Chemie GmbH
✓	Kohrsolin	BODE Chemie GmbH
✓	Kohrsolin Extra 4%	BODE Chemie GmbH
✓	Cosimed 5%	cosiMed GmbH
✓	Desomed rapid	Desomed Dr. Trippen GmbH
✓	Cleanisept Wipes	Dr. Schumacher GmbH
✓	Perasafe	DuPont/ Berner Pro
✓	Virkon	DuPont/ Berner Pro
✓	FD366 (+FD360)	Dürr Dental AG
✓	Brial	Ecolab
✓	Incidin Plus 3%	Ecolab
✓	Klercide 60/40	Ecolab
✓	Incidin OxyWipe	Ecolab
✓	DesiPlint	Frescon Oy
✓	Clinell Spray	GAMA Healthcare
✓	Chlor-Clean	Guest Medical
✓	Antibac Overflatedesinfeksjon 75%	KiiltoClean AS (Antibac AS)
✓	Easydes	KiiltoClean Oy
✓	Klorilli	KiiltoClean Oy
✓	Umonium ³⁸	Laboratoire Huckert's International
✓	Microzid	Schülke & Mayr GmbH
✓	Wet wipe	Wet Wipe A/S
+	Aldasan	Lysoform Dr. hans Rosemann
+	Antifect plus	Schülke & Mayr GmbH
o	Surfanios Premium	Laboratoire Anios
o	Surfanios Premium	Laboratoire Anios

- ✓ = recommended
 - +
 - o
- + = recommended to use diluted, 0.50% dilution
 o = recommended to use diluted, 0.25% dilution